

## Anti-Mouse Agarose VHH Beads

### Product Information

- Catalog Number: KTSM1341
- Volume: 1.0 mL (50% anti-mouse IgG nanobody coated agarose beads)
- Storage condition: 4°C for 12 months.
- Binding capacity: 1.0 mL Anti-Mouse Agarose VHH Beads slurry binds to ~4 mg mouse IgG
- Species Reactivity: Mouse
- Reactivity: bind to mouse IgGs

### Product Description

Anti-Mouse Agarose VHH Beads are suspension of activated agarose beads coupled with anti-mouse IgG nanobodies. This product is suitable for pulling down protein that is coupled with mouse IgGs or has mouse Fc tag.

### Application

Immunoprecipitation (IP), co-immunoprecipitation (CoIP), chromatin immunoprecipitation (CHIP), RNA-binding protein immunoprecipitation (RIP), enzyme activity detection, mass spectrometry analysis, etc.

### Application Note

Before usage, please invert the vial several times (DO NOT VORTEX) to form beads suspension, and take out beads' suspension using a pipette tip with bigger opening (i.e., cut the tip of a 1000 µl pipette tip with sterile scissors). Once opened, it is recommended to use paraffin film to seal the cap.

### Other Information

- Beads size: ~45 µm (7.5% cross-linked agarose beads)
- Storage buffer: 1XPBS, 0.03% sodium azide, 50% glycerol.

## Working Procedures:

### Cell Lysate Preparation

- Add 50  $\mu$ l Anti-Mouse Agarose VHH Beads to a microcentrifuge tube loaded with 500  $\mu$ l cell lysate sample, and incubate on ice for 30 min.
- Spin at 4°C and 1,200 g for 3 min and transfer the supernatant to a new microcentrifuge tube.

### Immunoprecipitation

- Add 5  $\mu$ g primary antibody to the microcentrifuge tube containing the precleared lysate, incubate for 1 hr.
- Add 50  $\mu$ l Anti-Mouse Agarose VHH Beads. Incubate for 1 hr on a rocking platform shaker.
- Centrifuge the tube at 4°C and 1,200g for 3 min.
- Remove supernatant completely and wash the pelleted beads 3 times with 500  $\mu$ l Lysis Buffer (50 mM Tris HCl, pH 8.0; 150 mM NaCl; 1% NP-40).

### SDS-PAGE

- Remove the supernatant from the last washing
- Add 100  $\mu$ l Laemmli Buffer (with 50 mM DTT or 2%  $\beta$ -mercaptoethanol) to bead pellet.
- Vortex and boil at 90-100°C for 10 min.
- Spin at 1,200g for 3 min,
- Collect the supernatant,
- Load samples for SDS-PAGE analysis (devoid of any beads).

## Results

Mouse IgGs IP from mouse serum using Anti-Mouse Agarose VHH Beads, as shown in the SDS-PAGE results:



M: Marker; I: Input; FT: Flow through; W: Wash; E: Elution

**Disclaimer: Products are for life science research only. Not for use in diagnostic procedures unless otherwise indicated.**